


Anti-Fat10/Ubd Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15205

Item Number: CM0027287



Introduction

Rabbit polyclonal antibody targeting human FAT10/UBD (ubiquitin D). Suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human and mouse. Protein weight 18kDa. Ideal for studying ubiquitin-proteasome pathway, immunity, and renal diseases.

[Learn More](#)

Parameter	Specification
Product Name	FAT10/UBD Rabbit pAb
Remarks/Alias	FAT10; UBD-3; GABBR1; FAT10/UBD
Species	Human
Gene ID (Human)	10537
Gene ID	10537
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human FAT10/UBD (NP_006389.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15205
Protein Weight	18kDa
Shipping	Ice bag