

Anti-Fatp4/Slc27A4 Monoclonal Antibody For Wb, Elisa - Q6P1M0

Item Number: CM0007279

Introduction

Recombinant anti-FATP4/SLC27A4 rabbit monoclonal antibody validated for WB and ELISA, detecting human and mouse FATP4, a long-chain fatty acid transporter and acyl-CoA ligase involved in lipid metabolism. Suitable for research on fatty acid transport, metabolism, and skin barrier function.

[Learn More](#)

Parameter	Value
Product Name	FATP4/SLC27A4 Rabbit mAb
Remarks/Alias	IPS; FATP4; ACSVL4; FATP4/SLC27A4
Species	Human
GeneID (Human)	10999
GeneID	10999
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-237 of human FATP4/SLC27A4 (NP_005085.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6P1M0
Protein Weight	72kDa
Shipping	Ice bag