

Anti-Fatty Acid Synthase (Fasn) Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P49327

Item Number: CM0002607



Introduction

This rabbit monoclonal antibody targets human Fatty Acid Synthase (FASN), validated for WB, IHC-P, IF/ICC, IP, ELISA. Reacts with human, mouse, and rat FASN, suitable for fatty acid biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	Fatty Acid Synthase (FASN) Rabbit mAb
Remarks/Alias	FAS; OA-519; SDR27X1; Fatty Acid Synthase (FASN)
Species	Human
GeneID (Human)	2194
GeneID	2194
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 2400-2500 of human FASN (P49327).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49327
Protein Weight	273kDa
Shipping	Ice bag