

Anti-Fatty Acid Synthase (Fasn) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P49327

Item Number: CM0018101



Introduction

Rabbit polyclonal antibody against human FASN, validated for WB, IHC-P, IF/ICC, ELISA. Detects endogenous FASN in human, mouse, rat samples. Suitable for fatty acid metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Fatty Acid Synthase (FASN) Rabbit pAb
Remarks/Alias	FAS; OA-519; SDR27X1; Fatty Acid Synthase (FASN)
Species	Human
GeneID	2194
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 2221-2511 of human Fatty Acid Synthase (FASN) (NP_004095.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49327
Protein Weight	273kDa
Shipping	Ice bag