

Anti-Fbln5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ubx5

Item Number: CM0029158

Introduction

Rabbit polyclonal antibody against human Fibulin-5 (FBLN5) for elastic fiber and ECM research. Validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	FBLN5 Rabbit pAb
Remarks/Alias	EVEC; UP50; ADCL2; ARMD3; CMT1H; DANCE; ARCL1A; FIBL-5; HNARMD; FBLN5
Species	Human
GeneID (human)	10516
GeneID	10516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 199-448 of human FBLN5 (NP_006320.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBX5
Protein Weight	50kDa
Shipping	Ice bag