


Anti-Fbp1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P09467

Item Number: CM0006652



Introduction

High-specificity anti-FBP1 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Targets gluconeogenesis rate-limiting enzyme fructose-1,6-bisphosphatase 1 (P09467).

[Learn More](#)

Parameter	Value
Product Name	FBP1 Rabbit mAb
Remarks/Alias	FBP; FBP1
Species	Human
Gene ID (Human)	2203
Gene ID	2203
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FBP1 (P09467).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09467
Protein Weight	37kDa
Shipping	Ice bag