


Anti-Fbp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P09467

Item Number: CM0029328



Introduction

FBP1 Rabbit pAb (P09467) is a polyclonal antibody raised in rabbit against a synthetic peptide of human FBP1. Validated for WB, IF/ICC, ELISA, and cross-reacts with human, mouse, and rat. Ideal for gluconeogenesis and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	FBP1 Rabbit pAb
Remarks/Alias	FBP; FBP1
Species	Human
GeneID (Human)	2203
GeneID	2203
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human FBP1 (NP_000498.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09467
Protein Weight	37kDa
Shipping	Ice bag