


Anti-Fbp2 Polyclonal Antibody For Wb, Elisa - 000757

Item Number: CM0013855



Introduction

FBP2 (Fructose-1,6-bisphosphatase isozyme 2) Rabbit Polyclonal Antibody. Validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Useful for investigating glycogen synthesis from lactate in skeletal muscle and related metabolic pathways.

[Learn More](#)

Parameter	Value
Product Name	FBP2 Rabbit pAb
Remarks/Alias	CORLK; FBP2
Species	Human
Gene ID (Human)	8789
Gene ID	8789
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-339 of human FBP2 (NP_003828.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00757
Protein Weight	37kDa
Shipping	Ice bag