

Anti-Fbx32/Fbxo32 Polyclonal Antibody For Wb, If-P, Elisa - Q969P5

Item Number: CM0026320

Introduction

Rabbit polyclonal antibody against human Fbx32/FBXO32 (Atrogin-1, MAFbx), suitable for WB, IF-P, and ELISA. Recognizes human and mouse samples. A key tool for studying muscle atrophy and ubiquitin-mediated proteolysis.

[Learn More](#)

Parameter	Value
Product Name	Fbx32/FBXO32 Rabbit pAb
Remarks/Alias	Fbx32; MAFbx; Fbx32/FBXO32
Species	Human
Gene ID	114907
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 206-355 of human Fbx32/FBXO32 (NP_478136.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q969P5
Protein Weight	42kDa
Shipping	Ice bag