

Anti-Fbx32/Fbxo32 Rabbit Monoclonal Antibody For Wb, Elisa - Q969P5

Item Number: CM0005016

Introduction

High-quality rabbit monoclonal antibody against human Fbx32/FBXO32 (Atrogin-1/MAFbx), an E3 ubiquitin ligase involved in muscle atrophy. Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	Fbx32/FBXO32 Rabbit mAb
Remarks/Alias	Fbx32; MAFbx; Fbx32/FBOX32
Species	Human
Gene ID (Human)	114907
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human Fbx32/FBOX32 (Q969P5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q969P5
Protein Weight	42kDa
Shipping	Ice bag