

Anti-Fbxl12 Polyclonal Antibody For Wb And Elisa - Q9Nxk8

Item Number: CM0015872

Introduction

High-quality rabbit polyclonal antibody targeting human FBXL12 (F-box/LRR-repeat protein 12). Validated for WB and ELISA, cross-reacts with mouse. Essential tool for studying ubiquitin-mediated proteolysis and cell cycle arrest.

[Learn More](#)

Parameter	Value
Product Name	FBXL12 Rabbit pAb
Remarks/Alias	Fbl12; FBXL12
Species	Human
Gene ID (Human)	54850
Gene ID	54850
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-326 of human FBXL12 (NP_060173.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NXK8
Protein Weight	37kDa
Shipping	Ice bag