


Anti-Fbxl2 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Ukc9

Item Number: CM0016723



Introduction

Rabbit polyclonal antibody targeting human FBXL2 (FBL2/FBL3), suitable for Western blot and ELISA in human, mouse, and rat samples. Ideal for research on ubiquitin-mediated proteolysis, cell cycle regulation, and inflammatory signaling.

[Learn More](#)

Parameter	Value
Product Name	FBXL2 Rabbit pAb
Remarks/Alias	FBL2; FBL3; FBXL2
Species	Human
GeneID (Human)	25827
GeneID	25827
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-240 of human FBXL2 (NP_036289.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKC9
Protein Weight	47kDa
Shipping	Ice bag