

Anti-Fbxl5 Polyclonal Antibody For Wb, Elisa - Q9Uka1

Item Number: CM0029812



Introduction

High-quality rabbit polyclonal antibody targeting FBXL5 (F-box/LRR-repeat protein 5). Validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Suitable for research on iron homeostasis, ubiquitin-proteasome system, and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	FBXL5 Rabbit pAb
Remarks/Alias	FBL4; FBL5; FLR1; FBXL5
Species	Human
GeneID (Human)	26234
GeneID	26234
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-310 of human FBXL5 (NP_036293.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKA1
Protein Weight	79kDa
Shipping	Ice bag