

Anti-Fbxo10 Rabbit Polyclonal Antibody For Wb And Elisa - Q9UK96

Item Number: CM0020105

Introduction

This rabbit polyclonal antibody targets human FBXO10, a component of the SCF E3 ubiquitin ligase complex. Validated for WB and ELISA with cross-reactivity to rat. Suitable for apoptosis and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	FBXO10 Rabbit pAb
Remarks/Alias	FBX10; PRMT11; FBXO10
Species	Human
GeneID (human)	26267
GeneID	26267
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human FBXO10 (NP_036298.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9UK96
Protein Weight	105kDa
Shipping	Ice bag