

Anti-Fbxo2 Polyclonal Antibody For Wb And Elisa - Q9UK22

Item Number: CM0029037

Introduction

FBXO2 rabbit polyclonal antibody validated for WB and ELISA with mouse and rat cross-reactivity. Targets F-box protein 2, a substrate recognition component of the SCF ubiquitin ligase involved in glycoprotein ER-associated degradation.

[Learn More](#)

Parameter	Value
Product Name	FBXO2 Rabbit pAb
Remarks/Alias	FBG1; FBX2; Fbs1; OCP1; NFB42; FBXO2
Species	Human
Gene ID (Human)	26232
Gene ID	26232
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-296 of human FBXO2 (NP_036300.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UK22
Protein Weight	33kDa
Shipping	Ice bag