

Anti-Fbxo21 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O94952

Item Number: CM0020006

Introduction

High-quality anti-FBXO21 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat FBXO21. Ideal for ubiquitination and protein degradation research.

[Learn More](#)

Parameter	Value
Product Name	FBXO21 Rabbit pAb
Remarks/Alias	FBX21; FBXO21
Species	Human
Gene ID (Human)	23014
Gene ID	23014
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 440-621 of human FBXO21 (NP_055817.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O94952
Protein Weight	72kDa
Shipping	Ice bag