


Anti-Fbxo22 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nez5

Item Number: CM0025459



Introduction

High-quality rabbit polyclonal antibody specific for human FBXO22 (FBX22, FISTC1). Validated for Western blot and ELISA applications; detects a 45kDa band. Ideal for ubiquitin pathway and muscle biology research.

[Learn More](#)

Parameter	Value
Product Name	FBXO22 Rabbit pAb
Remarks/Alias	FBX22; FISTC1; FBXO22
Species	Human
GeneID (human)	26263
GeneID	26263
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 104-403 of human FBXO22 (NP_671717.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NEZ5
Protein Weight	45kDa
Shipping	Ice bag