

Anti-Fbxo28 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nvf7

Item Number: CM0022449

Introduction

Rabbit polyclonal antibody targeting human FBXO28 (UniProt Q9NVF7), validated for WB and ELISA. Reacts with human and rat FBXO28. Suitable for ubiquitin ligase, genomic stability, and beta-cell survival studies.

[Learn More](#)

Parameter	Value
Product Name	FBXO28 Rabbit pAb
Remarks/Alias	Fbx28; DEE100; CENP-30; FBXO28
Species	Human
GenelD (Human)	23219
GenelD	23219
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NVF7
Protein Weight	41kDa
Shipping	Ice bag