

Anti-Fbxo3 Polyclonal Antibody For Wb, Elisa - Q9Uk99

Item Number: CM0017621



Introduction

Rabbit polyclonal antibody targeting human FBXO3 (F-box only protein 3), validated for Western blot (WB) and ELISA. Shows cross-reactivity with mouse. Suitable for studying E3 ubiquitin ligase function and protein degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	FBXO3 Rabbit pAb
Remarks/Alias	FBA; FBX3; FBXO3
Species	Human
Gene ID (Human)	26273
Gene ID	26273
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human FBXO3 (NP_036307.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UK99
Protein Weight	55kDa
Shipping	Ice bag