

Anti-Fbxo31 Polyclonal Antibody For Wb - Q5Xux0

Item Number: CM0012914



Introduction

Rabbit polyclonal antibody targeting human FBXO31 (aa 240-539), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying DNA damage and cell cycle regulation. Protein weight ~61 kDa.

[Learn More](#)

Parameter	Value
Product Name	FBXO31 Rabbit pAb
Remarks/Alias	FBX14; Fbx31; MRT45; FBXO14; pp2386; FBXO31
Species	Human
GeneID (Human)	79791
GeneID	79791
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 240-539 of human FBXO31 (NP_079011.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q5XUX0
Protein Weight	61kDa
Shipping	Ice bag