

Anti-Fbxo32 (Fbx32, Mafbx) Rabbit Polyclonal Antibody - Q969P5

Item Number: CM0013424

Introduction

Rabbit polyclonal antibody against human FBXO32 (Fbx32, MAFbx). Validated for IF/ICC and ELISA. Cross-reacts with mouse and rat. Suitable for skeletal muscle atrophy and ubiquitination studies.

[Learn More](#)

Parameter	Value
Product Name	Fbx32/FBXO32 Rabbit pAb
Remarks/Alias	Fbx32; MAFbx; Fbx32/FBXO32
Species	Human
Gene ID (Human)	114907
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 206-355 of human Fbx32/FBXO32 (NP_478136.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	Q969P5
Protein Weight	42kDa
Shipping	Ice bag