

Anti-Fbxo32 Rabbit Polyclonal Antibody For If/Icc, If-P, Elisa - Q969P5

Item Number: CM0026561

Introduction

High-quality rabbit polyclonal antibody against human FBXO32 (Atrogin-1, MAFbx). Validated for IF/ICC, IF-P, ELISA; reacts with human, mouse, rat. Ideal for muscle atrophy and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	FBXO32 Rabbit pAb
Remarks/Alias	Fbx32; MAFbx
Species	Human
GeneID (Human)	114907
GeneID	114907
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 201-300 of human FBXO32 (NP_478136.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969P5
Protein Weight	42kDa
Shipping	Ice bag