

Anti-Fbxo4 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukt5

Item Number: CM0028171



Introduction

FBXO4 Rabbit Polyclonal Antibody suitable for WB and ELISA, detects human and mouse FBXO4. The antibody targets F-box only protein 4, a substrate recognition subunit of SCF E3 ubiquitin ligase, useful for studying ubiquitin-mediated degradation.

[Learn More](#)

Parameter	Value
Product Name	FBXO4 Rabbit pAb
Remarks/Alias	FBX4; FBXO4
Species	Human
GeneID (Human)	26272
GeneID	26272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 138-387 of human FBXO4 (NP_036308.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UKT5
Protein Weight	44kDa
Shipping	Ice bag