

Anti-Fbxo42 Polyclonal Antibody For Wb - Q6P3S6

Item Number: CM0013173

Introduction

Rabbit polyclonal antibody targeting FBXO42 (JFK, Fbx42). Validated for WB, IHC-P, ELISA in human, mouse, rat. Useful for studying SCF ubiquitin ligase and p53 degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	FBXO42 Rabbit pAb
Remarks/Alias	JFK; Fbx42; FBXO42
Species	Human
Gene ID (Human)	54455
Gene ID	54455
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 468-717 of human FBXO42 (NP_061867.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6P3S6
Protein Weight	78kDa
Shipping	Ice bag