

Anti-Fbxw11 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukb1

Item Number: CM0026827



Introduction

High-quality rabbit polyclonal antibody against human FBXW11 (beta-TrCP2), validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat. Suitable for ubiquitin-proteasome, Wnt, and NF- κ B signaling studies.

[Learn More](#)

Parameter	Value
Product Name	FBXW11 Rabbit pAb
Remarks/Alias	Hos; BTRC2; FBW1B; Fbw11; BTRCP2; FBXW1B; NEDJED; FBXW11
Species	Human
Gene ID (human)	23291
Gene ID	23291
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 443-542 of human FBXW11 (NP_036432.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKB1
Protein Weight	62kDa
Shipping	Ice bag