

## Anti-Fbxw11 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Ukb1

Item Number: CM0021466



### Introduction

FBXW11 Rabbit pAb is a rabbit polyclonal antibody against human FBXW11. Suitable for WB, IF/ICC, ELISA and reacts with human, mouse, rat. Ideal for ubiquitination, Wnt signaling, and NF-κB pathway studies.

[Learn More](#)

| Parameter        | Value                                                    |
|------------------|----------------------------------------------------------|
| Product Name     | FBXW11 Rabbit pAb                                        |
| Remarks/Alias    | Hos; BTRC2; FBW1B; Fbw11; BTRCP2; FBXW1B; NEDJED; FBXW11 |
| Species          | Human                                                    |
| Gene ID (Human)  | 23291                                                    |
| Immunogen        | A synthetic Peptide of human FBXW11.                     |
| Source           | Rabbit                                                   |
| Category         | Polyclonal Antibodies                                    |
| Application      | WB, IF/ICC, ELISA                                        |
| Cross Reactivity | Human, Mouse, Rat                                        |
| SWISS            | Q9UKB1                                                   |
| Protein Weight   | 62kDa                                                    |
| Shipping         | Ice bag                                                  |