

Anti-Fbxw2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ukt8

Item Number: CM0024498

Introduction

Anti-FBXW2 rabbit polyclonal antibody validated for Western blotting and ELISA. Detects human, mouse, and rat FBXW2, a substrate-recognition component of the SCF E3 ubiquitin ligase complex. Suitable for ubiquitin pathway studies.

[Learn More](#)

Parameter	Value
Product Name	FBXW2 Rabbit pAb
Remarks/Alias	Md6; FBW2; Fwd2; FBXW2
Species	Human
Gene ID (Human)	26190
Gene ID	26190
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 102-454 of human FBXW2 (NP_036296.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UKT8
Protein Weight	52kDa
Shipping	Ice bag