

Anti-Fbxw5 Polyclonal Antibody For Wb, Elisa - Q969U6

Item Number: CM0028987

Introduction

Human FBXW5 rabbit polyclonal antibody validated for Western blot and ELISA. Targets the F-box/WD repeat-containing protein 5 (64 kDa), a substrate recognition component of SCF and DCX E3 ubiquitin ligase complexes—ideal for ubiquitination, cell cycle, and centriole duplication research.

[Learn More](#)

Parameter	Value
Product Name	FBXW5 Rabbit pAb
Remarks/Alias	Fbw5; FBXW5
Species	Human
Gene ID (Human)	54461
Gene ID	54461
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 184-566 of human FBXW5 (NP_061871.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q969U6
Protein Weight	64kDa
Shipping	Ice bag