

Anti-Fbxw7 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q969H0

Item Number: CM0030890



Introduction

Rabbit polyclonal antibody against FBXW7, suitable for Western blot, IHC (paraffin), and ELISA. Detects ~80 kDa FBXW7 in human, mouse, and rat samples. Targets the ubiquitin ligase substrate recognition component for proteasomal degradation studies.

[Learn More](#)

Parameter	Value
Product Name	FBXW7 Rabbit pAb
Remarks/Alias	AGO; CDC4; FBW6; FBW7; hAgo; FBX30; FBXW6; SEL10; hCdc4; DEDHIL; FBXO30; SEL-10; FBXW7
Species	Human
GeneID (human)	55294
GeneID	55294
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 231-371 of human human (NP_361014.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q969H0
Protein Weight	80kDa
Shipping	Ice bag