

Anti-Fbxw8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N3Y1

Item Number: CM0022138



Introduction

Rabbit polyclonal antibody targeting FBXW8, a substrate-recognition component of the Cul7-RING ubiquitin ligase complex. Validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Useful for ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	FBXW8 Rabbit pAb
Remarks/Alias	FBW6; FBW8; FBX29; FBXW6; FBXO29; FBXW8
Species	Human
GeneID (Human)	26259
GeneID	26259
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 500-598 of human FBXW8 (NP_699179.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N3Y1
Protein Weight	67kDa
Shipping	Ice bag