


Anti-Fcar Polyclonal Antibody For Wb, If/Icc, Elisa - P24071

Item Number: CM0027373



Introduction

Rabbit polyclonal antibody against human FCAR (CD89), validated for WB, IF/ICC, ELISA.

Recognizes human and mouse, ideal for studying Fc receptor function in myeloid cells and mucosal immunity.

[Learn More](#)

Parameter	Value
Product Name	FCAR Rabbit pAb
Remarks/Alias	CD89; FcalphaR; FcalphaRI; CTB-61M7.2; FCAR
Species	Human
GeneID (human)	2204
GeneID	2204
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-287 of human FCAR (NP_001991.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24071
Protein Weight	32kDa
Shipping	Ice bag