


Anti-Fcer1 Alpha Monoclonal Antibody For Wb, Elisa - P12319

Item Number: CM0007649



Introduction

High-affinity monoclonal antibody against human FcεR1 alpha (FCER1A), validated for WB and ELISA. Suitable for detecting human, mouse, and rat IgE receptor in allergy and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	FcεR1 alpha Rabbit mAb
Remarks/Alias	FCε1A, FcεRI, FCER1A, FcεR1 alpha
Species	Human
Gene ID (Human)	2205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-205 of human FcεR1 alpha (NP_001992.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12319
Protein Weight	30 kDa
Shipping	Ice bag