


Anti-FcεR1A Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P12319

Item Number: CM0012678



Introduction

High-affinity IgE receptor subunit alpha (FcεR1α) rabbit monoclonal antibody validated for WB, flow cytometry, and ELISA; reacts with human; ideal for allergy and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	FcεR1α Rabbit mAb
Remarks/Alias	FCE1A; FcERI; FCERIA
Species	Human
GeneID (Human)	2205
GeneID	2205
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 26-205 of human FcεR1α (NP_001992.1).)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human
SWISS	P12319
Protein Weight	30kDa
Shipping	Ice bag