


Anti-Fcf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y324

Item Number: CM0023145



Introduction

Rabbit polyclonal against human FCF1, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human and rat. Targets ribosome biogenesis, localized to nucleolus.

[Learn More](#)

Parameter	Value
Product Name	FCF1 Rabbit pAb
Remarks/Alias	Bka; UTP24; CGI-35; C14orf111; FCF1
Species	Human
Gene ID	51077
Immunogen	Synthetic peptide. A synthetic peptide corresponding to a sequence within amino acids 100-198 of human FCF1 (NP_057046.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9Y324
Protein Weight	23kDa
Shipping	Ice bag