

Anti-Fcgr1B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q92637

Item Number: CM0021046

Introduction

Rabbit polyclonal antibody against human FCGR1B (CD64b), validated for WB, IHC-P, and ELISA. Targets the putative high affinity IgG Fc receptor involved in humoral immunity.

[Learn More](#)

Parameter	Value
Product Name	FCGR1B Rabbit pAb
Remarks/Alias	CD34; FCG1; FcRI; CD64b; FCGR1; IGFR1; IGFRB; FCGR1B; FcgammaRIa
Species	Human
GeneID (human)	2210
GeneID	2210
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human FCGR1B (NP_001017986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q92637
Protein Weight	32kDa
Shipping	Ice bag