


Anti-Fcgr2A Polyclonal Antibody For Wb, If-P, Elisa - P12318

Item Number: CM0015369



Introduction

Rabbit polyclonal antibody targeting human FCGR2A (CD32), validated for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse and rat. Suitable for studying Fc receptor-mediated phagocytosis and immune responses.

[Learn More](#)

Parameter	Value
Product Name	FCGR2A Rabbit pAb
Remarks/Alias	CD32; FCG2; FcGR; CD32A; CDw32; FCGR2; IGFR2; FCGR2A1; FcgammaRIIa; FCGR2A
Species	Human
GeneID (Human)	2212
GeneID	2212
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 34-316 of human FCGR2A (NP_001129691.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12318
Protein Weight	35kDa
Shipping	Ice bag