

Anti-Fcgr3A Polyclonal Antibody For Wb, Elisa - P08637

Item Number: CM0024900



Introduction

Rabbit polyclonal antibody against human FCGR3A (CD16), validated for WB and ELISA, cross-reactive with mouse and rat. Suitable for studying antibody-dependent cellular cytotoxicity and NK cell biology.

[Learn More](#)

Parameter	Value
Product Name	FCGR3A Rabbit pAb
Remarks/Alias	CD16, FCG3, CD16A, FCGR3, IGFR3, IMD20, FCR-10, FCRIII, CD16-II, FCGRIII, FCRIIIA, FcGRIIIA, FCGR3A
Species	Human
GeneID (human)	2214
GeneID	2214
Immunogen	Synthetic peptide corresponding to amino acids 150-250 of human FCGR3A (NP_000560.7)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08637
Protein Weight	29kDa
Shipping	Ice bag