

Anti-Fcgr3B Rabbit Polyclonal Antibody For Wb And Elisa - O75015

Item Number: CM0026711



Introduction

High-quality rabbit polyclonal antibody targeting human FCGR3B (Fc gamma receptor III-B). Validated for Western blot (WB) and ELISA applications, with cross-reactivity to mouse and rat. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Specification
Product Name	FCGR3B Rabbit pAb
Remarks/Alias	CD16; FCG3; CD16A; CD16b; FCGR3; CD16-I; FCGR3A; FCR-10; FCRIII; FCRIIIb; FCGR3B
Species	Human
GeneID (Human)	2215
GeneID	2215
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-233 of human FCGR3B (NP_000561.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75015
Protein Weight	26kDa
Shipping	Ice bag