

Anti-Fcgrt/Fcrn Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P55899

Item Number: CM0007531



Introduction

Rabbit monoclonal antibody targeting FCGR1/FCRN (FcRn), suitable for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Recombinant immunogen. Ideal for investigating IgG transport and recycling.

[Learn More](#)

Parameter	Value
Product Name	FCGR1/FCRN Rabbit mAb
Remarks/Alias	FCRN; FcγRn; IgG Fc fragment receptor transporter alpha chain; FCGR1/FCRN
Species	Human
Gene ID (Human)	2217
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-290 human FCGR1/FCRN (NP_004098.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55899
Protein Weight	40kDa
Shipping	Ice bag