

Anti-Fcho2 Polyclonal Antibody For Wb And Elisa - Q0Jrz9

Item Number: CM0017494



Introduction

Rabbit polyclonal antibody against human FCHO2, validated for Western blot and ELISA. Targets an early endocytosis protein involved in clathrin-coated pit formation; cross-reacts with mouse. Ideal for studying membrane trafficking.

[Learn More](#)

Parameter	Value
Product Name	FCHO2 Rabbit pAb
Remarks/Alias	FCHO2
Species	Human
Gene ID (Human)	115548
Gene ID	115548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human FCHO2 (NP_620137.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q0JRZ9
Protein Weight	89kDa
Shipping	Ice bag