


Anti-Fcn1 Rabbit Polyclonal Antibody For Wb And Elisa - 000602

Item Number: CM0030182



Introduction

Rabbit polyclonal antibody against human FCN1 (M-ficolin). Suitable for Western blot and ELISA. Recognizes human, mouse, and rat. Ideal for innate immunity and complement pathway research.

[Learn More](#)

Parameter	Value
Product Name	FCN1 Rabbit pAb
Remarks/Alias	FCNM; FCN1
Species	Human
GeneID (Human)	2219
GeneID	2219
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-326 of human FCN1 (NP_001994.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	000602
Protein Weight	35kDa
Shipping	Ice bag