

Anti-Fcn2 Polyclonal Antibody For Wb, Elisa - Q15485

Item Number: CM0019187



Introduction

Rabbit polyclonal antibody against human FCN2 (Ficolin-2), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for innate immunity and complement pathway studies.

[Learn More](#)

Parameter	Value
Product Name	FCN2 Rabbit pAb
Remarks/Alias	P35; FCNL; EBP-37; ficolin-2; FCN2
Species	Human
GeneID (human)	2220
GeneID	2220
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 230-290 of human FCN2 (NP_004099.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15485
Protein Weight	34kDa
Shipping	Ice bag