

Anti-Fcrl3 Rabbit Polyclonal Antibody For Wb And Elisa - Q96P31

Item Number: CM0014351

Introduction

Rabbit polyclonal antibody against human FCRL3 (MAIA), validated for WB and ELISA. Suitable for research on B-cell biology, immune regulation, and fertilization. Reacts with human and mouse FCRL3.

[Learn More](#)

Parameter	Value
Product Name	FCRL3 Rabbit pAb
Remarks/Alias	MAIA; FCRH3; IFGP3; IRTA3; SPAP2; CD307c; FCRL3
Species	Human
Gene ID (Human)	115352
Gene ID	115352
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 18-180 of human FCRL3 (NP_443171.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96P31
Protein Weight	81kDa
Shipping	Ice bag