


Anti-Fcrl5 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Rd9

Item Number: CM0014698



Introduction

High-quality rabbit polyclonal antibody targeting human FCRL5 (FcRH5, CD307), validated for WB and ELISA. Ideal for B-cell development, lymphoma research, and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	FCRL5 Rabbit pAb
Remarks/Alias	CD307; FCRH5; IRTA2; BXMAS1; CD307e; PRO820; FCRL5
Species	Human
Gene ID (Human)	83416
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 16-190 of human FCRL5 (NP_001182317.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96RD9
Protein Weight	106kDa
Shipping	Ice bag