

Anti-Fdps Rabbit Polyclonal Antibody For Wb, Elisa - P14324

Item Number: CM0027968



Introduction

Rabbit polyclonal antibody against FDPS (FPPS), validated for WB and ELISA. Reacts with human, mouse, rat. Targets farnesyl pyrophosphate synthase, a key enzyme in isoprenoid biosynthesis. Suitable for cancer, metabolic, and lipid signaling research.

[Learn More](#)

Parameter	Value
Product Name	FPS/FDPS Rabbit pAb
Remarks/Alias	FPS; FPPS; POROK9; FPS/FDPS
Species	Human
Gene ID (Human)	2224
Gene ID	2224
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 120-419 of human FPS/FDPS (NP_001995.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14324
Protein Weight	48kDa
Shipping	Ice bag