

Anti-Fdx1 Polyclonal Antibody For Wb And Elisa - P10109

Item Number: CM0024183



Introduction

This rabbit polyclonal antibody targets human FDX1 (Adrenodoxin), a mitochondrial ferredoxin essential for steroidogenesis. Validated for WB and ELISA, it cross-reacts with mouse and is ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	FDX1 Rabbit pAb
Remarks/Alias	FDX1; ADX; FDX; LOH11CR1D; ferredoxin 1
Species	Human
Gene ID (Human)	2230
Gene ID	2230
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 61-160 of human FDX1 (NP_004100.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10109
Protein Weight	19kDa
Shipping	Ice bag