

Anti-Fdx1/Adx Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P46656

Item Number: CM0029696

Introduction

Rabbit polyclonal antibody against ferredoxin 1 (FDX1/ADX), suitable for WB, IP, IF/ICC, IHC-P, and ELISA. Detects human, mouse, and rat. Ideal for steroidogenesis and mitochondrial electron transport research.

[Learn More](#)

Parameter	Value
Product Name	FDX1/ADX Rabbit pAb
Remarks/Alias	FDX1; ADX; FDX; LOH11CR1D; ferredoxin 1
Species	Mouse
GeneID (Human)	2230
GeneID	14148
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 99-188 of mouse FDX1/ADX (NP_032022.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46656
Protein Weight	20 kDa
Shipping	Ice bag