

Anti-Fdx1/Adx Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10109

Item Number: CM0003531



Introduction

Rabbit monoclonal antibody targeting human FDX1/ADX (P10109), essential for steroid hormone synthesis. Validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	FDX1/ADX Rabbit mAb
Remarks/Alias	ADX; FDX; LOH11CR1D; FDX1/ADX
Species	Human
GeneID (Human)	2230
GeneID	2230
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FDX1/ADX (P10109).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10109
Protein Weight	19 kDa
Shipping	Ice bag