

Anti-Fdx1/Adx Rabbit Polyclonal Antibody For Wb, Elisa - P10109

Item Number: CM0024736



Introduction

High-quality rabbit polyclonal antibody targeting human FDX1/ADX (adrenodoxin). Validated for WB and ELISA with human samples. Investigate steroidogenesis, mitochondrial electron transport, and adrenal gland biology. UniProt P10109.

[Learn More](#)

Parameter	Value
Product Name	FDX1/ADX Rabbit pAb
Remarks/Alias	ADX; FDX; LOH11CR1D; FDX1/ADX
Species	Human
GeneID (Human)	2230
GeneID	2230
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human FDX1/ADX (NP_004100.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10109
Protein Weight	19kDa
Shipping	Ice bag