

Anti-Fdxr Rabbit Polyclonal Antibody For Wb, Elisa - P22570

Item Number: CM0026273

Introduction

FDXR Rabbit Polyclonal Antibody validated for WB and ELISA, reactive with human, mouse, rat. Suitable for mitochondrial P450 system studies, coenzyme Q biosynthesis research, and redox metabolism investigations.

[Learn More](#)

Parameter	Value
Product Name	FDXR Rabbit pAb
Remarks/Alias	ADR; ADXR; ANOA; FDXR
Species	Human
GeneID (Human)	2232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 252-491 of human FDXR (NP_077728.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22570
Protein Weight	54kDa
Shipping	Ice bag